

CIO Applications europa

CIOAPPLICATIONSEUROPE.COM

Two new products for your Valab® solution.

By Eric Rogari, CEO, VALAB



Eric Rogari

(turnaround time). To date, our intelligent autoverification solution assists more of 360 structures across 11 countries in the hard decision-making process. To continue meeting the needs of laboratories, VALAB now offers from an additional Business Intelligence (BI) module (ValView), which complements its biological validation support solution, as well as Valab® license management tool (VSP).

ValView has been designed to simplify and enhance the monitoring of your laboratory's activity. This innovative solution offers a seamless and comprehensive user experience. With just a few clicks, you can get an overview of Valab® statistical analyses, thanks to several dashboards featuring dynamic graphs and interactive tables. Closer to Valab®'s day-to-day performance, the visualization of trends and variations in production indicators (Expertise rate, validation rate, etc.) over adjustable periods facilitates maintenance and ensures optimal ongoing control. Also available in ValView, the intuitive reading and easy interpretation of the evolution of various Valab® rates, accompanied

The review of results in the laboratory total testing process is a key step before the release of the medical reports. Valab®, a solution highly sensitive, robust, and sustainable, supports physicians on a daily basis, in avoiding particularly the release of false-negative. Indeed, Valab® ensures an efficient and secure autoverification of medical reports by filtering and alerting for critical and debatable values. This leads to harmonized medical validation, reduced errors, saved time and faster TAT



Valab aims to be an ally to medical laboratories in their validation process, which is why our Valab expert system and complementary products provide precise, reliable, and highly sensitive autoverification, while offering flexibility, customization, and security. Through increasingly advanced automation and standards, the VSP and ValView respectively ensure the continuity of services and strengthen and simplify quality monitoring for automated validation by Valab, always better aligned with the needs of medical laboratories today

by their detailed explanations (Rejection criteria), enable instant evaluation of the performance of the Valab® system.

In a more in-depth manner, ValView also allows you to see the ranking of the most frequently blocked analyses, as well as the drift in validation rates per analysis, with the option to adjust the selection parameters. At a glance, the key points that need optimization in Valab® can be identified, enabling the expert system to be calibrated according to the specific needs of the laboratory, thus enhancing the efficiency of the automatic validation performed by Valab®.

In short, managing your laboratory's activity with VALAB's new tool has never

been easier. ValView provides flexibility, real-time reactivity, ease of analysis and optimal security in monitoring your laboratory's activity, helping you to gain control and efficiency while complying with ISO 15189 standards.

The VSP (VALAB Services Provider) is a hosted license server that offers seamless, simplified license management packed with valuable benefits. The service continuity is ensured with the VSP, which also improves the responsiveness of VALAB teams via a monitoring and an alert system linked directly to our teams. Moreover, users have more flexibility with their Valab® license, as the number of files autoverified per day can be temporarily exceed up to 12 times per year.

These two products meet the growing need for more efficient quality control and monitoring, ensuring compliance with regulatory standards that guarantee optimized and safe laboratories processes. **CA**

Today, we can see that technology is transforming innovation at its core, allowing organizations to test new ideas at a faster pace and pocket-friendly prices—that were unimaginable even a decade ago. At an era where one can stick features on web sites and tell within hours how customers respond, decision makers need a platform that can tell them which technology can be beneficial for their businesses and increase their ROI. CIO Applications Europe has fast culminated into the most sought after magazine Pan-Europe as it has transformed its course hand-in-hand with the technology and has culminated into a leading media brand, with a unique editorial focus on bringing to light the core innovations in technology.

CIO Applications Europe magazine stands out with its unique approach of learning from peers approach offering professionals the most comprehensive collection of technology trends. The magazine is the single most successful initiative Pan-Europe to advice and guide decision makers regarding the latest in the fast evolving technology landscape. We are determined to propose a myriad of additional services that can improve businesses and help customers deal with issues related to this industry.