



System Requirements Valab® Client V14.xx



Characteristic	Possible values
Usage	Fat client for remote access to the Valab® server
OS	Windows 10® Windows Server 2016® x64 Windows Server 2012® R2 x64
Processor	dual core - 64 bits - 2 Ghz
RAM	4 Gbytes
Disk space	Space necessary for the OS + 20 Gbytes
Graphics card	1366 * 768 minimum recommended
.NET Framework	Version 4.5.1 minimum
Setup configuration / FireWall - Network	Open ports 2110, 8089, 8090 and 8080 to the Valab® server
Access	Read / write access to the installation directory of the Valab® client for the user