



VALAB Services Provider
(Hosted License Server)

**For enhanced management
of your Valab[®] license,
the VSP improves**

01

Flexibility

**Authorized
overrun**

02

Reactivity

Monitoring

03

Security

**Hosting and
access**

04

Robustness

**Service
continuity**

www.valab.com

With VSP, a more flexible license, an even more responsive, robust service and secure hosting

01. A more flexible license

Thanks to VSP, the management of your Valab® license is optimized to maintain your laboratory's activities, without any service interruptions.



Authorized overrun

- **12 days** per year
- Up to **twice the number of medical files** per day



Remote administration

Installation, license modification, management of addendums or readjustment of your offer, etc.

VSP makes it easy and secure to manage your licenses remotely, without having to access your site via remote maintenance.

03. Secure hosting

The access and transfer of VSP data are fully secured, ensuring a highly reliable service.



Data encryption*:
exchanged data is encrypted via **SSL protocol** and **protected by token**.



Authenticated communication:
the **Valab® server** always **initiates the exchange** with the **VSP**.



Secure hosting:
the **VSP** is **hosted on a high-security server**.

* No confidential or personal data is collected by the VSP.

02. A more responsive service

VSP monitoring enhances the responsiveness of VALAB Support.



Alerts

- **Alerts sent to VALAB support:**

License level exceeded, no files appraised (SIL connection stopped), appraisal rate and/or file validation too low, end of license validity and/or temporary amendment, loss of connection between the VSP and your Valab®, etc.



Indicators

- **Statistical counters:** NRR, NER and NVR
- **Information :** installed version, automation, consultative mode, constant monitoring, etc.



04. A more robust service

To ensure your service is always up and running, and to manage your license quickly and simply.



Service continuity:

- **Robustness: no impact** following OS or server changes; compatible with server farms or high-availability clusters.
- **Stability:** in the event of a technical problem, your VSP license adapts to limit service interruptions, etc.

VSP connection

Connection status: Active Customer ID: [redacted] Refresh

Server status: Nominal Server ID: [redacted]

License

Selected licence: 53 - main - G40-BCGH-VAV-X60-C30-F30-(8) State: Assigned Expiry date: -

Limits

Reports/day: 4000 Users: 30

Licensed tests: 100 LIS connections: 1

Modules

BIOCHEMISTR COAGULATION

BLOOD GASES HAEMATOLOG

Options

ValView

No. of authorisations: 12 / 12

Resetting date: 8/29/2025

Active temporary riders

List of temporary riders active on the selected licence:

ID	Start date	End date	Reports/day	Tests	Users	LIS connection	Modules	Options
4	8/1/2024	1/1/2025	5000	-	-	-	-	-



Overview window:

View your license information at a glance!



Streamlined license management



Extended service continuity



Secure communications

Would you like more information on the benefits of VSP?

Don't hesitate to contact us!



contact@valab.com



+33 (0) 531 083 499



www.valab.com



01.09.2024